
50 10

1

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pH electrode 50 10

Introduction

The **50 10** is a *Low Cost* pH electrode for general use. Low maintenance electrode, as it does not require monitoring of electrolyte level. Suitable for routine measurements.

Electrode preparation

Remove the storage protector containing electrolyte by unscrewing it ③. Verify that there are not air bubbles inside the membrane ⑦. They can be removed by shaking the electrode as a clinical thermometer.

Calibration

⚠ See *pH-meter instruction manual*.

Measurements. Recommendations

- During calibration and measurement the buffers and the samples must be slightly stirred. The stirring speed should be identical in both cases.
- The sample must cover the electrode's diaphragm ⑥.
- After measurement, rinse the electrode with distilled water. Dry it with soft tissue, only by contact, without rubbing.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.

⚠ *Do not leave the electrode immersed in distilled water.*

Maintenance

- Maintain the membrane's surface clean ⑦. Protect it from friction or knocking.
- When an electrode gets contaminated, the simple rinsing with distilled water is not sufficient. The most suitable cleaning solution is one which acts in a more selective way over the contamination (look for HACH regeneration solutions).
- **Re-hydration.** A "dry" electrode recovers by immersing it in HACH cleaning electrode solution, code C20C370, or in diluted HCl solution for several hours. It is recommended to read the information supplied with the electrode's regenerating solutions.

Troubleshooting

... the pH-meter always measures the same pH value.

Short circuit in the electrode or connector.

Substitute the electrode.

... unstable reading.

Air bubble in the membrane.

Shake the electrode as a thermometer.

... impossible to perform calibration.

... very slow response.

Dirt on the diaphragm or the membrane.

De-hydrated membrane.

Electrode ageing or wear.

Clean, regenerate or substitute the electrode.

Storage

Place the electrode in the storage protector ③ with KCl 3M.

Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not respond appropriately after regeneration, must be replaced by a new one.

 *If an electrode has too short life span, less than three months, ask us about it. Probably there is a more suitable electrode for the application.*

Specifications

pH measuring range	0 ... 14
Operating temp.	0 ... 80 °C
Reference element	Ag wire coated with AgCl
Diaphragms	ceramic
Electrolyte	gel
Body material	external: polycarbonate, internal: glass
Min. immersion depth	12 mm 

Applications

Aqueous samples in general.

Limitations

Low conductivity, dirty or viscous solutions.

Products with colloids or solids in suspension.

Samples containing sulphides, reducing sugars or other substances which react with the silver ion.

Temperatures exceeding 80 °C.

Warranty

HACH guarantees this electrode against manufacturing defects.

Validity: 6 months from shipment date.

Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

QUALITY CERTIFICATE

The attached electrode 50 10 has passed the quality test and follows the next specifications:

- Asymmetry potential < ± 15 mV.
- Sensitivity, pH 4...7 (at 25 °C) > 98%.
- Response time, pH 4...7 < 20 s.

HACH COMPANY World Headquarters

P.O. Box 389, Loveland, CO 80539-0389 U.S.A.

Tel. (970) 669-3050

(800) 227-4224 (U.S.A. only)

Fax (970) 669-2932

orders@hach.com

www.hach.com

HACH LANGE GMBH

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D-40549 Düsseldorf, Germany

Tel. +49 (0) 2 11 52 88-320

Fax +49 (0) 2 11 52 88-210

info@de.hach.com

www.de.hach.com

HACH LANGE Sàrl

6, route de Compois

1222 Vérenaz

SWITZERLAND

Tel. +41 22 594 6400

Fax +41 22 594 6499