
52 08

pH

DOC022.98.90537
09-2015, Edition 1

Puntos clave. *Punti chiave.* Points clé. *Key points.*

① Tapón protector.
Tappo protettore.
Bouchon protecteur.
Protective cap.

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Introduction

The **52 08** is a pH electrode for microsamples. It can take measurements in volumes under 100 µl.

⚠ Warning. Extremely fragile electrode. Handle with care.

Electrode preparation

Remove the storage protector containing electrolyte (6). Verify that there are not air bubbles inside the membrane (8). They can be removed by shaking the electrode as a clinical thermometer. Remove the cap from the electrolyte filling port (3).

Calibration

⚠ See pH-meter instruction manual.

Measurements. Recommendations

- During calibration and measurement the buffers and the samples must be slightly stirred. The stirring speed should be identical in both cases.
- The sample must cover the electrode's diaphragm (7).
- After measurement, immediately rinse the electrode with distilled water. Dry it with soft tissue, only by contact, without rubbing.
- Between measurements the electrode's membrane should be immersed in aqueous solution. The most appropriate place is the storage protector with some electrolyte.

⚠ Do not leave the electrode immersed in distilled water.

Maintenance

- Periodically refill the electrode with electrolyte. Maintain the level as high as possible (4).
- Maintain the membrane's surface clean (8). Protect it from friction or knocking.
- When an electrode gets contaminated, the simple rinsing with distilled water is not sufficient. The most suitable cleaning solution is one which acts in a more selective way over the contamination (look for HACH regeneration solutions).
- **Re-hydration.** A "dry" electrode recovers by immersing it in HACH cleaning electrode solution or in diluted HCl solution for several hours. It is recommended to read the information supplied with the electrode's regenerating solutions.

Troubleshooting

... always measures the same pH value.
Short circuit in the electrode or connector.

Substitute the electrode.

... unstable reading.

Air bubble in the membrane.

Shake the electrode as a thermometer.

... impossible to perform calibration.

... very slow response.

Dirt on the diaphragm or the membrane.

De-hydrated membrane.

Electrode ageing or wear.

Clean, regenerate or substitute the electrode.

Storage

The connector should be protected from humidity by the protective cap ①.

Place the electrode in the storage protector ⑥ with KCl 3M.

Close the refilling port with the stopper ③.

Life of an electrode

The average life expectancy of a pH electrode is approximately one year depending on the working conditions.

An electrode used for more than six months that does not respond appropriately after regeneration, must be replaced by a new one.

 *If an electrode has too short lifespan, less than three months, ask us about it. Probably there is a more suitable electrode for the application.*

Specifications

| | |
|----------------------|--|
| Measuring range | 0... 14 |
| Operating temp. | 0... 80 °C |
| Reference element | encapsulated Ag/AgCl crystals |
| Diaphragms | ceramic |
| Electrolyte | KCl 3M |
| Body material | glass |
| Min. immersion depth | 5 mm  |
| Diameter | 3 mm |

Applications

Biology and clinical areas.

Limitations

Viscous solutions or samples with colloids.
Distilled water.

Warranty

HACH guarantees this electrode against manufacturing defects.

Validity: 6 months from shipment date.

Limitations:

- Accidental damage.
- Inadequate applications.
- Non-fulfilment of the recommendations.
- Normal wear and tear of the electrode.

QUALITY CERTIFICATE

The attached electrode 52 08 has passed the quality test and follows the next specifications:

- Asymmetry potential $< \pm 15$ mV.
- Sensitivity, pH 4...7 (at 25 °C) $> 98\%$.
- Response time, pH 4...7 < 20 s.

HACH COMPANY World Headquarters

P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
Tel. (970) 669-3050
(800) 227-4224 (U.S.A. only)
Fax (970) 669-2932
orders@hach.com
www.hach.com

HACH LANGE GMBH

Willstätterstraße 11
D-40549 Düsseldorf, Germany
Tel. +49 (0) 2 11 52 88-320
Fax +49 (0) 2 11 52 88-210
info@de.hach.com
www.de.hach.com

HACH LANGE Sàrl

6, route de Compois
1222 Vérenaz
SWITZERLAND
Tel. +41 22 594 6400
Fax +41 22 594 6499