



Be Right™



Pipet Filler, Safety Bulb

Réf. produit:	1465100
EUR Prix:	Nous contacter
Disponible	

Unique cone adapter allows universal pipet filling. Simply squeeze rubber bulb to charge for filling, place pipet into cone adapter and insert into sample. Release the bulb to fill. Remove bulb and cover pipet end with finger for dispensing. Easy and safe. Autoclavable.

Caractéristiques

Volume / taille du carton: 1 each

Contenu de la boîte

Includes autoclavable elastic wristband.