



Be Right™



## Buffer Solution, pH 7.00 (NIST), color-coded yellow, 4L

Réf. produit: 2283556

**ARTICLE OBSOLETE**  
Cet article n'est plus disponible.

### pH Buffer standard solutions for laboratory and field analyses

Solution tampon de qualité, code couleur jaune pour identification aisée. Traçabilité NIST. pH  $7,00 \pm 0,02$  à 25 °C.

For pH meter standardization. Color-coded for easy identification. pH is  $\pm 0.02$  at 25°C. 4 Liters

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Caractéristiques

Code couleur:	Oui
Concentration:	$7,00 \pm 0,02$ pH
Conditions de stockage:	10 - 25 °C (protect from light)
Description:	Buffer Solution (color-coded Yellow) (NIST)
Intervalle d'entretien:	24 months
Notes spéciales:	Observed shelf-life after opening, following recommended storage and use: 6 months
Paramètres:	pH
Propriétés du contenant :	Bottle
Traçable NIST:	Oui
Valeurs pH:	7.00
Volume / taille du carton:	4 l