



Be Right™



Bandelettes de test pour phosphate, 50 unités.

Réf. produit: 2757150
EUR Prix: Nous contacter
Disponible

Test strips in a variety of parameters.

High-quality Hach chemistry on an easy-to-use test strip. Simply dip the strip in sample and compare with colors printed on the bottle to find a reading. Test strips are accurate, portable, and low-cost. Applications include drinking water, wastewater, pool and spa, industrial process water, and more. Parameters available include chlorine, nitrate and nitrite, hardness, alkalinity, chloride, ammonia, and more. Hach offers a high-value 5 in 1 test strip, which measures free chlorine, total chlorine, hardness, alkalinity, and pH - all on one strip.

Facile à utiliser, jetable et très abordable

Un moyen idéal pour obtenir des résultats chiffrés et rapides sur le terrain ou en laboratoire

Caractéristiques

Modèle de kit de test:	Bandelettes de test
Paramètres:	Phosphore, orthophosphate (réactif) - Selon la teneur en PO_4
Plage de mesure:	0 - 50 ppm PO_4
Plateforme:	Test Strip
Plus petites incréments:	Incréments : 0/5/15/30/50
Poids de livraison:	0.2