



**Be Right™**



## **Pipet, serological, sterile, polystyrene w/ plug, 2 mL, 50/pk**

<b>Réf. produit:</b>	<b>2092826</b>
EUR Prix:	Nous contacter
Livraison en 2 semaines	

Sterile, disposable plastic pipet with plug. Bulk packaged in 25/pack Color coded. Black ascending/descending graduations. Calibrated "To Deliver". Accuracy  $\pm 2\%$ .

---

### **Caractéristiques**

Capacité:	2 mL
Graduations:	0.01 mL
Volume / taille du carton:	50 /pk