



Be Right™



BAG, WHIRL-PAK* 5-1/2X9 PK/10

Réf. produit:	1437297
EUR Prix:	Nous contacter
Livraison en 2 semaines	

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection.

Caractéristiques

Capacité:	720 mL
Volume / taille du carton:	10