



Be Right™



m-Endo broth, dehydrated, 500 g

Réf. produit:

1462300

Dangereux

EUR Prix:

Nous contacter

Livraison en 2 semaines



Hach microbiological chemicals

Hach carries a wide variety of media and accessories for microbiology applications, including dehydrated media, nutrient agar, broths, acids, dilution water, sterilization indicator, and reagents. Combined with Hach's microbiology expertise, these products mean making your own media in the lab for standard or custom sampling applications is economical and practical. Instructions for mixing are included on each bottle. Dehydrated media product is available for your application, including m-Endo, Nutrient, Tryptic, m-Green, and many others.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Caractéristiques

Volume / taille du carton:

500 g