



Be Right™



DEHA Color Disc Test Kit, Model DH-1

Réf. produit: 2168200

ARTICLE OBSOLETE
Cet article n'est plus disponible.

Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. The Model DH-1 test kit is a boiler/cooling water kit utilizing the color disc method. N, N-Diethylhydroxylamine (DEHA) replaces hydrazine and its derivatives as a oxygen scavenger in steam generation systems. Contains 50 tests.

Simple, fast, color disc method

Measures N, N-Diethylhydroxylamine (DEHA) directly in three different concentration ranges

To measure other oxygen scavenger compounds, use Prod. No. 2349500, Oxygen Scavenger Test Kit

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Caractéristiques

| | |
|------------------------|-------------------------------|
| Modèle de kit de test: | DH-1 |
| Nom de la méthode: | Color disc/ Iron reduction |
| Nombre de tests: | 100 |
| Paramètres: | DEHA |
| Plage de mesure: | 1 - 65 µg/L DEHA |
| Plateforme: | Disque de couleur |
| Poids de livraison: | 9 |
| Smallest increment: | 1 µg/L 5 µg/L 22.5 µg/L |
| Type de boîtier: | R |