



Glycol Test Kit (Model EG-1)

Réf. produit: 2186400
EUR Prix: Nous contacter
Livraison en 2 semaines

Semiquantitative determination of glycol in water or oil. Visual appearance of colour indicates glycol presence.

Ethylene glycol (anti-freeze) contamination of automotive lubricants is an indication of a cracked engine block. Glycol presence also can indicate leaks in cooling systems. Test Kit Model EG-1 presents a purple color if ethylene glycol and other 1,2 glycols are present in oil or water. The kit contains approximately 25 tests and comes with a 14.4" x 8" x 7" Blue Polypropylene Carrying Case.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Caractéristiques

Modèle de kit de test:	EG-1
Nom de la méthode:	Visual/ Appearance of color
Nombre de tests:	25
Paramètres:	Glycol - Ethylene glycol in oil or water
Plage de mesure:	Presence or absence of glycol, 0- greater than 150 mg/L
Poids de livraison:	3.5
Type de boîtier:	G