



## Test kit, glutaraldehyde, model GT-1

Réf. produit: 2587200  
EUR Prix: Nous contacter  
Demander un délai de livraison

From oil fields to cooling towers, hospitals to food plants, poultry farms to dental offices, the Glutaraldehyde Test Kit can help monitor disinfection processes. The test kit measures from 0.5 to 4000 ppm glutaraldehyde.

**Helps monitor disinfection processes**

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Caractéristiques

Alimentation:	NA
Échelle:	1 - 2800 mg/L multiple ranges
Modèle:	GT-1
Modèle de kit de test:	GT-1
Nom de la méthode:	Color disc
Nombre de tests:	100
Paramètres:	Glutaraldehyde
Plage de mesure:	NA
Plateforme:	Disque de couleur
Poids de livraison:	8.3
Smallest increment:	0.5
Type de boîtier:	Q