



## Hardness (Total & Calcium) Test Kit, Model HA-4P MG-L

Réf. produit: 145701  
EUR Prix: Nous contacter  
Disponible

Hardness, Total & Calcium (Magnesium by difference) drop count titration/EDTA test kit, range 20 - 400 mg/L CaCO<sub>3</sub>, 100 tests.

The Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field.

**Utilises drop count titration - a simple, economical titration method**

---

### Caractéristiques

Modèle de kit de test:	HA-4P MG-L
Nom de la méthode:	Drop count titration/EDTA
Nombre de tests:	100
Paramètres:	Hardness, total & Calcium - (Magnesium by difference)
Plage de mesure:	20 - 400 mg/L CaCO <sub>3</sub>
Plateforme:	Titrateur goutte à goutte
Smallest increment:	20
Type de boîtier:	D

---

### Contenu de la boîte

Includes five reagents, bottle, measuring tube, instruction sheet, and carrying case.