



Be Right™



Hardness (Total & Calcium) Test Kit, Model HAC-DT

Réf. produit: 2063900
EUR Prix: Nous contacter
Demander un délai de livraison

Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10 - 4000 mg/L. Contains approximately 100 of each (Total and Calcium) test.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Caractéristiques

Modèle de kit de test:	HAC-DT
Nom de la méthode:	Digital Titrator/EDTA
Nombre de tests:	100 Ca +100 T
Paramètres:	Hardness, Total & Calcium - as CaCO ₃
Plage de mesure:	10 - 4,000 mg/L CaCO ₃
Plateforme:	Titrateur digital
Poids de livraison:	4
Smallest increment:	0.1-10
Type de boîtier:	F

Contenu de la boîte

Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.