



## Free & Total Chlorine Test Kit (Model CN-70)

Réf. produit:

1454200

Dangereux

EUR Prix:

Nous contacter

Disponible



Simple, low-cost color disc method is the most accurate visual test kit method

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Caractéristiques

Modèle de kit de test:	CN-70 <sup>2</sup>
Nom de la méthode:	Color disc/DPD
Nombre de tests:	100 Free, 100 Total
Note de pied de page:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Paramètres:	Chlorine, low range - As free & total Cl <sub>2</sub>
Plage de mesure:	0 - 3,4 mg/L Cl <sub>2</sub>
Plateforme:	Disque de couleur
Poids de livraison:	2
Smallest increment:	0.02
	0.1
Type de boîtier:	D