



Be Right™



Dissolved Oxygen AccuVac® Color Disc Test Kit, 0-1,000 µg/L

Réf. produit: 2501050

ARTICLE OBSOLETE
Cet article n'est plus disponible.

Remplacements suggérés:

- 2501025

Low DO levels usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This kit uses the fast, accurate visual method for determining dissolved oxygen concentrations. Range (µg/L): 0-1,000. Kit contains 25 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

AccuVac® vials mean no mess, no measuring

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Caractéristiques

Modèle de kit de test:	AccuVac® Kit
Nom de la méthode:	Color disc/LRDO
Nombre de tests:	25
Paramètres:	Dissolved Oxygen - As DO
Plage de mesure:	0 - 1,000 µg/L O ₂
Poids de livraison:	4.1
Smallest increment:	50 µ g/L
Type de boîtier:	G