



## Test kit, sulfate, model SF-1

Réf. produit:

225100

Dangereux

EUR Prix:

Nous contacter

Livraison en 2 semaines



Widely found in natural waters, sulfate may be at high levels in mine drainage. High concentrations of magnesium or sodium sulfate in water act as cathartics. This test kit provides economical determination of sulfate levels using the Extinction/Turbidimetric methods.

### Economical test kit for determination of sulfate levels using the Extinction/Turbidimetric Method

---

#### Caractéristiques

Modèle de kit de test:	SF-1
Nom de la méthode:	Extinction/Turbidimetric
Nombre de tests:	100
Paramètres:	Sulfate - As $\text{SO}_4^{2-}$
Plage de mesure:	50 - 200 mg/L
Poids de livraison:	1.5
Smallest increment:	50
Type de boîtier:	D