



Be Right™



Alkalinity test kit AL-AP, 5-100/20-400 mg/L CaCO₃ 100 pcs

Réf. produit: 2444301
EUR Prix: Nous contacter
Disponible

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists.

Method: Drop count titration

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Kit utilizes drop count titration method and contains 100 tests.

Simple, economical titration method

Rugged kit contains all required reagents, solutions, and glassware

Caractéristiques

Modèle de kit de test:	AL-AP, MG-L
Nom de la méthode:	Drop count titration/Sulfuric acid
Nombre de tests:	100
Note de pied de page:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Paramètres:	Alkalinity - (P) & Total (MO) as CaCO ₃
Plage de mesure:	5 - 100 mg/L
Plateforme:	Titrateur goutte à goutte
Poids de livraison:	1
Smallest increment:	5
	20
Type de boîtier:	B