



Be Right™



Nitrogen, Nitrite Colour Disc Test Kit, Model NI-6

Réf. produit:

224000

Dangereux

EUR Prix:

Nous contacter

Livraison en 2 semaines



Utilising the simple colour disc method/ferrous sulfate, this kit measures Nitrogen, Nitrite as NO_2^- and contains 100 tests. Range (mg/L): 0-8/0-1600.

Simple, low-cost colour disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Caractéristiques

Modèle de kit de test:	NI-6
Nom de la méthode:	Colour disc/ Ferrous sulphate
Nombre de tests:	100
Paramètres:	Nitrogen, Nitrite, high range - as NO_2^-
Plage de mesure:	0 - 8/0-1600 mg/L NO_2^-
Plateforme:	Disque de couleur
Poids de livraison:	2
Smallest increment:	2 40
Type de boîtier:	D

Contenu de la boîte

Includes colour disc, comparator box and tubes, reagents, deionised dilution water, instructions and carrying case.