



Be Right™



CEL/890 Advanced portable laboratory

Réf. produit:

2994000

Dangereux

EUR Prix:

Nous contacter

Demander un délai de livraison



Offering on-site results in an easy-to-use, rugged, portable lab, Hach's Oil & Gas Water Quality Laboratory covers parameters critical to effective source water and base water analysis. Depending on water composition, this kit may also be used to support analysis of produced water, flowback water, water treatment, and enhanced oil recovery applications.

Caractéristiques

Application:

Fracing

Contenu de la boîte

This laboratory includes: DR 900 Colorimeter, HQ40d Multi-Parameter Digital Meter, Digital Titrator, all required reagents for 12 critical parameters, apparatus and accessories, and two convenient, rugged carrying cases.