



Be Right™



Beaker, polypropylene, low form, 50 mL

Réf. produit:	108041
EUR Prix:	Nous contacter
Livraison en 2 semaines	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Caractéristiques

Capacité:	50 mL
Diamètre, extérieur:	51.1 mm with spout
Matériau:	Polypropylene
Matériel:	Polypropylene
Unités:	each