



**Be Right™**



## **Beaker, polypropylene, low form, 250 mL**

<b>Réf. produit:</b>	<b>108046</b>
EUR Prix:	Nous contacter
Disponible	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes  $\geq 250$  mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### **Caractéristiques**

Capacité:	250 mL
Diamètre, extérieur:	83.6 mm with spout
Matériau:	Polypropylene
Matériel:	Polypropylene
Unités:	each