



Be Right™



Test kit, RapidSilver™ visual

Réf. produit: 2674500
EUR Prix: Nous contacter
Demander un délai de livraison

Test Kit Silver, Model RapidSilver™. Range 0-50 ppb. Silver As Ag Visual Comparison.

Provides fast, accurate information on silver discharges

Used by surface finishers, photographic film manufacturers and processors, and others to prevent pollution and economic loss

Kit directions also include digestion and dilution procedures needed for testing silver strike bath samples

Caractéristiques

Alimentation:	NA
Échelle:	5 - 50 ppb Ag
Modèle:	NA
Modèle de kit de test:	RapidSilver™
Nom de la méthode:	Visual Comparison
Nombre de tests:	100
Paramètres:	Silver-As Ag
Plage de mesure:	NA
Poids de livraison:	4
Smallest increment:	5 µg/L
Type de boîtier:	F