



Be Right™



Cylinder, graduated, 100 mL +/- 0.6 mL, 1.0 mL divisions (white markings)

Réf. produit: 50842
EUR Prix: Nous contacter
Livraison en 2 semaines

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes ≥ 25 mL.

Caractéristiques

Capacité:	100 mL
Certifications:	Class B
Matériau:	Glass
Méthode de mesure:	To Deliver
Type:	Mesure unique
Volume / taille du carton:	1 each