



Be Right™



**Sension+ 5201 pH liquid filled combination electrode with automatic temperature compensation via built-in Pt1000 for general use. The encapsulated reference ensures high stability and long lifetime.**

Réf. produit: LZW5201.97.0002

**ARTICLE OBSOLETE**  
Cet article n'est plus disponible.

Remplacements suggérés:  
• LZW5014T.97.002

pH liquid electrode (refillable) with ceramic diaphragm and built-in A.T.C. Pt1000. Glass body. Range 0-14 pH, 0 to 100 °C. BNC and banana (Pt1000).

This pH electrode is suitable for aqueous media in general. The encapsulated reference ensures high stability and long lifetime.

### High stability and long lifetime

Encapsulated Ag/AgCl reference system

---

## Caractéristiques

Connecteurs:	BNC for pH + Banana for Pt 1000
Dimensions (D x L):	12 mm x 130 mm
Exactitude:	0.002 pH with pH31 Meter
Jonction de sonde:	Ceramic
Longueur:	130 mm
Matériel:	Sensor Body: Glass
Paramètres:	pH
Plage de mesure:	0 - 14 pH
Plage de température:	0 à 100 °C
Solution de remplissage:	LZW9500.99 KCl 3M
Type:	Lab
Type de capteur:	Glass
Type de sonde:	Laboratory